



INSTITUTE  
FOR RESEARCH  
IN BIOMEDICINE

# The Bio-PanPipe Software Package

---

Daniel Ortiz

February 13, 2019

# Table of Contents

---

1. Introduction
2. Package Overview
3. Main Tools and File Formats
4. Whole Pipeline Example

# Introduction

---

- Execution of bam file pipelines entails many difficulties:
  - Downloading of very large files
  - Combination of tools with different input requirements
  - Existence of dependencies between tools
  - Tools may need to be added or removed
  - Each tool has specific computational requirements
  - Pipeline may need to be executed for hundreds of files
  - Parallelism should be exploited when possible
  - ...
- Bio-PanPipe provides a PanPipe module as well as a set of utilities to tackle these problems

# Package Overview

---

# Package Dependencies

- Shell Bash
- Python
- Conda
- PanPipe
- Database download clients
- Slurm Workload Manager (optional)

# Package Installation

- Obtain the package using git:

```
git clone gitlab@fsupeksvr.irbbarcelona.pcb.ub.es:dortiz/bio-panpipe.git
```

- Change to the directory with the package's source code and type:

```
./reconf  
./configure  
make  
make install
```

**NOTE:** use `--prefix` option of `configure` to install the package in a custom directory

## Additional configure Options

- `--with-panpipe=DIR`: sets location PanPipe of package
- `--with-icgcstor=DIR`: sets location of ICGC's storage client
- `--with-aspera=DIR`: enables Aspera Connect download client
- `--with-egadecrypt=DIR`: location of EGA decryptor tool



- Execution of pipelines processing normal-tumor bam files
- Automate processing of all of the samples of a dataset
- Handle file downloading as part of pipeline execution

# Supported Databases and Download Clients

- Databases

- EGA
- ICGC

- Download clients

- aspc
- icgc-storage-client
  - Amazon cloud
  - Collaboratory cloud
- pyega3

# Implemented Analysis Steps

- **bam file downloading:**
  - `download_ega_{norm|tum}_bam`
  - `download_ega_asp_{norm|tum}_bam`
  - `download_aws_{norm|tum}_bam`
  - `download_collab_{norm|tum}_bam`
- **bam file manipulation:**
  - `sort_{norm|tum}_bam`
  - `index_{norm|tum}_bam`
  - `sambamba_mpileup_{norm|tum}_bam`
  - `delete_bam_files`

- **Small Indels and Single Nucleotide Variant Callers:**
  - `manta_germline`
  - `manta_somatic`
  - `platypus_germline`
  - `strelka_germline`
  - `strelka_somatic`

- Copy Number Variant Callers:

- `ascatngs`
- `cnvkit`
- `facets`
- `sequenza`
- `wisecondorx`

- **Structural Variant Callers:**
  - `delly`, `parallel_delly`
  - `lumpy`, `parallel_lumpy`
  - `parallel_svtyper`
- **MSI Analyzers:**
  - `msisensor`

# Main Tools and File Formats

---

- `query_ega_metadata`
- `query_icgc_metadata`
- `analyze_dataset`



- Extracts information from EGA metadata
- Main input parameters:
  - -s <string>: file with sample information
  - -a <string>: file with analysis information
  - -t <string>: file with study information
  - -p <string>: file listing Aspera box content
  - -f <int>: output format

- Extracts information from ICGC metadata
- Main input parameters:
  - -d <string>: file with donor information
  - -a <string>: file with aws manifest
  - -t <string>: table file in json format
  - -f <int>: output format:

- Uses metadata information to automate analysis of a whole dataset
- Main input parameters:
  - -pfile <string>: file with pipeline steps to be performed
  - -r <string>: file with reference genome
  - -m <string>: file with metadata, one entry per line

# The `bam_analysis.sh` Module

- Implements a PanPipe module for analyzing bam files
- Functions can be classified in 3 groups:
  - Download of bam files
  - Manipulation of bam files
  - Bioinformatics analysis (SNV, CNV and SV callers, MSI analyzers)

- Reference genome operations:
  - `filter_contig_from_genref`
  - `gen_bed_for_genome`
- Data preparation for analysis steps:
  - `convert_snppos_to_snpgcc`
  - `create_snv_pos_ascat`
  - `gen_wisecondorx_ref`

- **EGA/ICGC metadata:** information regarding a whole dataset that is typically spread out in a set of files
- **Analysis metadata:** file providing all the information of a given dataset that is relevant to automate its analysis
- **Analysis automation script:** file with a sequence of commands automating the analysis of a dataset

- Sample information (`Sample_File.map`)
  - contains file name info
- Analysis information (`Analisis_Sample_meta_info.map`)
  - contains donor and phenotype information
- Study information (`Study_analysis_sample.map`)
  - contains EGA id information
- Aspera box content (`dbbox_content`)

- Donor information (`donor.<study_name>.tsv`)
  - contains gender information
- AWS manifest (`manifest.aws-virginia.<code>.tsv`)
  - contains object id, file name and donor id information
- JSON table file (`icgc_table.json`)
  - contains phenotype information



# Analysis Metadata (EGA)

- Created with the `query_ega_metadata` tool
- Example entries:

```
EGAF00001664282 phenotype=Blood|Normal_blood gender=male ; EGAF00001664327 phenotype=Skin|  
  Tumour_metastasis_to_local_lymph_node gender=male  
  
EGAF00001670586 phenotype=Blood|Normal_blood gender=male ; EGAF00001664289 phenotype=Skin|  
  Tumour_metastasis_to_local_lymph_node gender=male  
  
EGAF00001664356 phenotype=Skin|Tumour_metastasis_to_distant_location gender=male ; EGAF00001670533  
  phenotype=Blood|Normal_blood gender=male  
  
EGAF00001661882 phenotype=Blood|Normal_blood gender=male ; EGAF00001661538 phenotype=Skin|  
  Tumour_metastasis_to_local_lymph_node gender=male  
...
```

# Analysis Metadata (EGA Aspera)

- Created with the query\_ega\_metadata tool
- Example entries:

```
EGAD00001003388/PART_2/EGAZ00001300436_20170516_AWS_MELA_3c3ed66c-1505-4614-ac9d-575a6713b06a.bam.crypt
  phenotype=Blood|Normal_blood gender=male ; EGAD00001003388/PART_3/
EGAZ00001300354_20170516_AWS_MELA_daf1ffd8-0a0f-4869-abc8-5be0b4fc1a21.bam.crypt phenotype=Skin|
Tumour_metastasis_to_local_lymph_node gender=male

EGAD00001003388/PART_3/EGAZ00001303407_20170516_AWS_MELA_a197619e-f3e2-41f6-aef7-d1fadb3c1f5b.bam.crypt
  phenotype=Blood|Normal_blood gender=male ; EGAD00001003388/PART_2/
EGAZ00001300389_20170516_AWS_MELA_3a9bf676-1a7b-4718-8396-fb36cc89b688.bam.crypt phenotype=Skin|
Tumour_metastasis_to_local_lymph_node gender=male

EGAD00001003388/PART_3/EGAZ00001300416_20170516_AWS_MELA_f64eba46-d8a1-46f2-ba66-1b509e16c946.bam.crypt
  phenotype=Skin|Tumour_metastasis_to_distant_location gender=male ; EGAD00001003388/PART_3/
EGAZ00001303394_20170516_AWS_MELA_7bb66858-7533-4f96-9cd4-41aae2fe18b2.bam.crypt phenotype=Blood|
Normal_blood gender=male

...
```

# Analysis Metadata (ICGC)

- Created with the `query_icgc_metadata` tool
- Example entries:

```
34fa2369-424f-5886-9d23-6d19f8f15278 tumor female ; d759d07f-330c-5d0c-bd28-af72147dfb17 normal female
284f1424-d250-59cf-b105-da277b061e4a normal female ; e7e69d23-fb0d-5d3d-9027-ebf355053dbf tumor female
c42fffad-4ffd-59ba-93f1-2c573547369c normal female ; 3a33ef20-dfd0-50b0-afc2-38de9a5baa32 tumor female
37f076d6-fa64-5b5d-a0d0-b5cd7428d4a2 normal female ; 2c34270b-98d2-54b9-bdd3-068c6a9d858f tumor female
...
```

# Pipeline Automation Script

- Created with the `analyze_dataset` tool (`-p` option)
- At each entry (one per line), PanPipe's `pipe_exec` tool is used to analyze a normal-tumor bam file pair
- Entry example:

```
/home/dortiz/bio/software/bam-utils/bin/pipe_exec --pfile /home/dortiz/bio/software/bam-utils/share/bam
-utils/examples/basic_test.ppl --outdir /mnt/raid/dortiz/bio/tasks/bam_analysis_testing_pipeline/
d759d07f-330c-5d0c-bd28-af72147dfb17_34fa2369-424f-5886-9d23-6d19f8f15278 --sched SLURM -r /home/
dortiz/bio/data/genome_references/refseq_hg19_filt.fa -extn d759d07f-330c-5d0c-bd28-af72147dfb17
-extt 34fa2369-424f-5886-9d23-6d19f8f15278 -g XX -cr /home/dortiz/bio/data/genome_references/
refseq_hg19_filt.fa.bed -sv /home/dortiz/bio/data/facets_info/00-common_all.vcf -sg /home/dortiz/
bio/data/ascatngs_info/r93/SnpGcCorrections_GRCh37_1000g.tsv -mc chrY -egastr 50 -egacred /home/
dortiz/bio/software/ega-download-client-python/dortiz_cred.json
```

# Extending Package Functionality

- Focus on `bam_analysis.sh` module
- Two mechanisms:
  - Add new functions directly in `bam_analysis.sh`
  - Define a complementary module and import it in addition to `bam_analysis.sh`

# Whole Pipeline Example

---

# Analysis File

```
#import bam_analysis
#
download_ega_norm_bam cpus=1 mem=2048 time=10:00:00 stepdeps=none
download_ega_tum_bam cpus=1 mem=2048 time=10:00:00 stepdeps=none
sort_norm_bam cpus=1 mem=4G time=10:00:00 stepdeps=afterok:download_ega_norm_bam
sort_tum_bam cpus=1 mem=4G time=10:00:00 stepdeps=afterok:download_ega_tum_bam
index_norm_bam cpus=1 mem=1G time=4:00:00 stepdeps=afterok:sort_norm_bam
index_tum_bam cpus=1 mem=1G time=4:00:00 stepdeps=afterok:sort_tum_bam
manta_somatic cpus=8 mem=3G time=6:00:00 stepdeps=afterok:index_norm_bam,afterok:index_tum_bam
strelka_somatic cpus=8 mem=6G time=6:00:00 stepdeps=afterok:index_norm_bam,afterok:index_tum_bam,
    afterok:manta_somatic
msisensor cpus=8 mem=6G time=5:00:00 stepdeps=afterok:index_norm_bam,afterok:index_tum_bam
facets cpus=1 mem=20G time=4:00:00 stepdeps=afterok:index_norm_bam,afterok:index_tum_bam
cnvkit cpus=8 mem=8G time=10:00:00 stepdeps=afterok:index_norm_bam,afterok:index_tum_bam
ascatngs cpus=8 mem=25G time=12:00:00 stepdeps=afterok:index_norm_bam,afterok:index_tum_bam
platypus_germline cpus=1 mem=4G time=5:00:00 stepdeps=afterok:index_norm_bam
delete_bam_files cpus=1 mem=1G time=0:10:00 stepdeps=afterok:manta_somatic,afterok:strelka_somatic,
    afterok:msisensor,afterok:cnvkit,afterok:facets,afterok:ascatngs,afterok:platypus_germline
```

# Pipeline

